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1. Das IJ, Cheng CW, Watts RJ, et al. Accelerator beam data commissioning equipment and procedures: report of the TG-106 of the Therapy Physics Committee of the AAPM. Med Phys. 2008 Sep;35(9):4186-215.
2. Chang Z, Wu Q, Adamson J, et al. Commissioning and dosimetric characteristics of TrueBeam system: composite data of three TrueBeam machines. Med Phys. 2012 Nov;39(11):6981-7018.
3. J Jang, YN Kang, HJ Shin, et al. Measurement of Beam Data for Small Radiosurgical Fields: Comparison of Cyberknife Multi-sites in Korea. Progr Nucl Sci Technol, 2011, 1: 537-540.
4. De Ost B, Schaeken B, Vynckier S, et al. Reference dosimetry for helical tomotherapy: practical implementation and a multicenter validation. Med Phys. 2011 Nov;38(11):6020-6.
5. Beyer GP. Commissioning measurements for photon beam data on three TrueBeam linear accelerators, and comparison with Trilogy and Clinac 2100 linear accelerators. J Appl Clin Med Phys. 2013 Jan 7;14(1):4077.
6. Glide-Hurst C, Bellon M, Foster R, et al. Commissioning of the Varian TrueBeam linear accelerator: a multi-institutional study. Med Phys. 2013 Mar;40(3):031719.
7. Kerns JR, Followill DS, Lowenstein J, et al. Technical Report: Reference photon dosimetry data

- for Varian accelerators based on IROC-Houston site visit data. *Med Phys.* 2016 May;43(5):2374.
8. Kerns JR, Followill DS, Lowenstein J, et al. Reference dosimetry data and modeling challenges for Elekta accelerators based on IROC-Houston site visit data. *Med Phys.* 2018 May;45(5):2337-2344.
  9. Akino Y, Mizuno H, Tanaka Y, et al. Inter-institutional variability of small-field-dosimetry beams among HD120™ multileaf collimators: a multi-institutional analysis. *Phys Med Biol.* 2018 Oct 18;63(20):205018.
  10. Tanaka Y, Mizuno H, Akino Y, et al. Do the representative beam data for TrueBeam™ linear accelerators represent average data? *J Appl Clin Med Phys.* 2019 Feb;20(2):51-62
  11. Das IJ, Ding GX, Ahnesjö A. Small fields: nonequilibrium radiation dosimetry. *Med Phys.* 2008 Jan;35(1):206-15.
  12. Das IJ, Njeh CF, Orton CG. Point/counterpoint: vendor provided machine data should never be used as a substitute for fully commissioning a linear accelerator. *Med Phys.* 2012 Feb;39(2):569-72.
  13. 日本医学物理学会編 外部放射線治療における水吸収線量の標準計測法 12.
  14. Klein EE, Hanley J, Bayouth, J et al. Task Group 142 report: quality assurance of medical accelerators. *Med Phys.* 2009 Sep;36(9):4197-212.
  15. Wen N, Zhao B, Kim J, et al. IMRT and RapidArc commissioning of a TrueBeam linear accelerator using TG-119 protocol cases. *J Appl Clin Med Phys.* 2014 Sep 8;15(5):4843.
  16. Szpala S, Cao F, Kohli K. On using the dosimetric leaf gap to model the rounded leaf ends in VMAT/RapidArc plans. *J Appl Clin Med Phys.* 2014 Mar 6;15(2):4484.
  17. Isono M, Akino Y, Mizuno H, et al. Inter-unit variability of multi-leaf collimator parameters for IMRT and VMAT treatment planning: a multi-institutional survey. *J Radiat Res.* In press.
  18. Ben M, Agnieszka O, Claudio F, et al. ESTRO BOOKLET NO. 7: Quality assurance of treatment planning systems. Practical examples for non-IMRT photon beams.
  19. Fraass B, Doppke K, Hunt M, et al. American Association of Physicists in Medicine Radiation Therapy Committee Task Group 53: quality assurance for clinical radiotherapy treatment planning. *Med Phys.* 1998 Oct;25(10):1773-829.
  20. 日本医学物理学会 タスクグループ 01 X 線治療計画システムに関する QA ガイドライン  
Quality Assurance of Treatment Planning Systems for X-ray Beams.

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